OM-C0353E 006

Contra Angle Hanpiece

Ti-Max X

OPERATION MANUAL

X10L / X10 / X12L / X12 / X15L / X15 / X20L / X20
X25L / X25 / X85L / X85 / X95L / X95

0197

MADE IN JAPAN

1. User and Indications for Use

User: Qualified Professionals

Indications for Use: This medical device is intended for transmitting rotation of the power source at the direct drive ratio or at a different gear ratio, thereby running the instrument such as a bur or a reamer to cut and polish natural or artificial teeth during dental treatment.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure. (Refer to "7. Post-use Maintenance")
- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components and may lead to blood coagulation inside the handpiece. Coagulated blood may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, NSK recommends to immediately lubricate those handpieces with PANA SPRAY Plus / PANA SPRAY.
- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue.
- Supply coolant water and coolant air while using the handpiece. No supplying the coolant water and coolant air may lead to overheating, causing burn injuries or product failure.
- Do not use the handpiece for prophylaxis purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or cause burn injuries due to the overheating of the handpiece head.
- Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.
- Ball bearings, in time, will become worn and will require replacement. Be sure to inspect the handpiece before use (Refer to "6. Checking the Handpiece Before Each Use"). If abnormal vibration or noise are found, the bearings or other internal parts may require replacement. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- Read this Operation Manual before use to fully understand the product functions and file for future reference.
- When operating the product always consider patient safety.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to "6. Checking the Handpiece Before Each Use")
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur.
 - Bur which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur.
- Do not use burs other than specified in "9. Specifications" as such burs may be accidentally released or break during rotation.
- Always keep the bur shank clean. Debris in the chuck could cause poor bur concentricity or low chuck retention force.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the handpiece head bearings, or accidental release of the bur, may occur.
- Do not exceed Max. Bur Length shown in "9. Specifications".
- Always follow the instructions provided by the bur manufacturer.

CAUTION

- Do not exceed the bur speed recommended by the bur manufacturer.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- The handpiece is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks. (Refer to "8. Periodical Maintenance Checks")
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

X85L / X85

- Use a short shank FG bur for the product because the handpiece has the miniature head. Using a long bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur.

X95L / X95

- Do not use a surgical bur with its maximum length of 25mm or longer. Using such bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur.

3. Accessory List

No.	Part Name	Quantity
1	E-Type Spray Nozzle	1

4. Connecting and Disconnecting the Handpiece from the Motor

4-1 Connecting

- 1) Insert the handpiece direct into the motor (Fig. 1). When connecting an optic handpiece, turn the handpiece until it clicks and locks into position.
- 2) Confirm that the handpiece is securely connected to the motor.

4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.


Fig. 1

⚠ CAUTION

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964 (EN ISO 3964)).
- Do not exceed the Max. Rotation Speed (Motor) shown in "9. Specifications".

5. Inserting and Removing the Bur

5-1 Inserting the Bur

X10L / X10 / X12L / X12 / X15L / X15 / X25L / X25 (Fig. 2)

- 1) Insert the bur until it is correctly secured in place.
- 2) Depress the Push Button and insert the bur into the chuck until the bur "notch" mechanism engages. Release the button.
- 3) Confirm that the bur is secure by gently pulling and pushing the bur without depressing the Push Button.

X20L / X20 / X85L / X85 / X95L / X95 (Fig. 3)

- 1) Insert the bur into the chuck.
- 2) Depress the Push Button to open the chuck (❶).
- 3) Insert the bur fully into the chuck until it stops (❷) then release the Push Button.
- 4) Confirm that the bur is secure by gently pulling and pushing the bur without depressing the Push Button.


Fig. 2


Fig. 3

NOTICE

- For easy operation of the Chuck Release Button simply grip the handpiece and use your thumb tip to press the Chuck Release Button.

5-2 Removing the Bur

Depress the Chuck Release Button to open the chuck (❶) and, while the chuck is open, remove the bur (❷).


Fig. 4

⚠ CAUTION

- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur.
 - Bur which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur.

CAUTION

- Do not use burs other than specified in "9. Specifications" as such burs may be accidentally released or break during rotation.
- Always keep the bur shank clean. Debris in the chuck could cause poor bur concentricity or low chuck retention force.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the handpiece head bearings, or accidental release of the bur, may occur.
- Do not exceed Max. Bur Length shown in "9. Specifications".
- Always follow the instructions provided by the bur manufacturer.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Do not insert or remove the bur until rotation has completely stopped.
- Do not apply excess pressure to the bur as it may break or bend and become difficult to remove.

6. Checking the Handpiece Before Each Use

Follow the check procedure below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1) Check the Head Cap is firmly tightened.
- 2) Check coolant water is flowing properly.
- 3) Insert the bur. (Refer to "5. Inserting and Removing the Bur")
- 4) Rotate the handpiece for about one minute, using coolant water, and with the motor operating at full speed. During rotation, check for abnormalities such as abnormal rotation, vibration, noise.
- 5) After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally (Fig. 5).


Fig. 5

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

7. Post-use Maintenance

After each patient maintain the product as follows.

WARNING

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure.

CAUTION

- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.

7-1 Preparation

- 1) Wear eye protection, a mask and gloves to avoid infection.
- 2) Remove the bur.
- 3) Disconnect the handpiece from the motor.

7-2 Cleaning

■ Manual Cleaning (External Surface)

- 1) Remove debris on the handpiece surface using a brush under running water (at 38°C or less).
DO NOT use a wire brush.
- 2) Wipe off moisture on the surface using an absorbent cloth.

CAUTION

- If water enters the handpiece, using an air syringe, blow air into the rear of the handpiece to remove the internal water.

■ Automatic Cleaning (External and Internal Surfaces)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).
Refer to the thermo-disinfector manual.

CAUTION

- After washing with thermo-disinfector, prior to lubrication, dry the product until all internal moisture is completely removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

7-3 Cleaning of Optic Illumination Points

Wipe clean all Optic Illumination Points (Glass Rod) using an alcohol-immersed cotton swab. Remove all debris and oil.


Fig. 6

CAUTION

- Do not use a sharp tool to clean the Glass Rod. It could damage the glass and reduce the light transmission.

7-4 Disinfection

■ Manual Disinfection (External Surface)

Wipe the surface with cloth soaked in disinfectant solution such as alcohol.

■ Automatic Disinfection (External and Internal Surfaces)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883). Refer to the thermo-disinfector manual.

7-5 Lubrication

■ NSK PANA SPRAY Plus / PANA SPRAY

Apply NSK PANA SPRAY Plus / PANA SPRAY every time after each use and/or before autoclaving.

- 1) Attach the CORRECT Spray Nozzle onto the spray port nozzle on the can.
- 2) Insert the Spray Nozzle into the rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds (Fig. 7).


Fig. 7

⚠ CAUTION

- Hold the spray can upright.
- Firmly hold the handpiece to prevent it slipping when spray pressure is applied.
- Spray lubricant until it expels from the handpiece head.

⚠ CAUTION

- If you wish to purge excessive oil from inside the handpiece, rotate the handpiece for approx 15 seconds without a bur. During the rotation, do not depress the chuck release button, especially when using absorbent cloth to prevent oil from scattering. Depressing the push button during rotation may damage the chuck mechanism.


Fig. 8


Fig. 9

NOTICE

- NSK recommends the use of "Spray Mist Absorber" (Y900084) to prevent oil mist expelling out of the handpiece head.

■ Chuck cleaning (X20L / X20 / X85L / X85 / X95L / X95)

Clean the Push Button chuck once a week.

- 1) Mount the Tip Nozzle onto the spray can port.
- 2) Lubricate the chuck directly where the bur is inserted (Fig. 10).
- 3) Lubricate handpiece by using NSK PANA SPRAY Plus / PANA SPRAY (Fig. 7) or NSK automatic handpiece cleaning and lubrication system.


Fig. 10

! CAUTION

- If the chuck is not regularly cleaned the chuck grip may be weakened and the bur may be accidentally released during use.

■ NSK automatic handpiece cleaning and lubrication system

When using NSK automatic handpiece cleaning and lubrication system, refer to the system instructions.

WARNING

- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components and may lead to blood coagulation inside the handpiece. Coagulated blood may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, NSK recommends to immediately lubricate those handpieces with PANA SPRAY Plus / PANA SPRAY.

7-6 Sterilization

Sterilize the product by autoclave sterilization. After each patient sterilize as noted below.

- 1) Insert into an autoclave pouch. Seal the pouch.
- 2) Autoclavable under the conditions below.
Autoclave for more than 20min. at 121°C, or 15min. at 132°C, or 3min. at 134°C.
- 3) The handpiece should remain in the autoclave pouch until required for use.

CAUTION

- Clean and lubricate the handpiece prior to autoclaving. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- To avoid product failure, do not use a sterilizer that exceeds a cycle temperature of 138°C, including the dry cycle. In some sterilizers, the chamber temperature may exceed 138°C. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.

NOTICE

- NSK recommends Class B sterilizers as stated in EN 13060.

8. Periodical Maintenance Checks

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Head cap is loose	Check that the Head Cap is firmly tightened.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.
Coolant Water	Operate the handpiece and check that the coolant water is flowing through all spray ports.

9. Specifications

Model	X10L	X10	X12L	X12	X15L	X15	X20L	X20	X25L	X25	X85L	X85	X95L	X95
Max. Rotation Speed (Motor)	40,000min ⁻¹													
Max. Rotation Speed (Handpiece)	2,500min ⁻¹		4,000min ⁻¹		10,000min ⁻¹		40,000min ⁻¹				200,000min ⁻¹			
Gear Ratio	16:1 Reduction		10:1 Reduction		4:1 Reduction		1:1 Direct Drive				1:5 Increasing			
Bur Type	ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm CA Bur						ISO 1797-1 (EN ISO 1797-1) Type3 Ø1.59-1.60mm Standard FG Bur		ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm CA Bur		ISO 1797-1 (EN ISO 1797-1) Type3 Ø1.59-1.60mm Short Shank FG Bur		ISO 1797-1 (EN ISO 1797-1) Type3 Ø1.59-1.60mm Standard FG Bur	
Chuck Length	11.6mm						11.9mm		11.6mm		10.2mm		11.6mm	
Max. Bur Length	22.5mm						25mm (Recommended 19mm)		22.5mm		21mm		25mm (Recommended 19mm)	
Max. Working Part Diameter	Ø4.0mm						Ø2.0mm		Ø4.0mm		Ø2.0mm			
Optic	Glass Rod	-	Glass Rod	-	Glass Rod	-	Glass Rod	-	Glass Rod	-	Glass Rod	-	Glass Rod	-
Water Spray Type	Single												Quattro	
Water Pressure	0.08 - 0.20MPa (0.8 - 2.0kgf/cm ²)													
Chip Air Pressure	0.10 - 0.20MPa (1.0 - 2.0kgf/cm ²)													

Water Consumption	Min. 50mL/min (0.2MPa)
Chip Air Consumption	Min. 1.5L/min (0.2MPa)
Use Environment	Temperature: 10 - 40°C, Humidity: 30 - 75% (No Condensation)
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa

* Available models differ according to countries/regions.

10. Symbol

This product is Autoclavable up to Max.135°C.

This product can be washed via Thermo Disinfector.

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Manufacturer.

Authorized representative in the European community.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

11. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

12. Spare Parts List

Model	Order Code
E-Type Spray Nozzle	Z019090

13. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK Europe GmbH