**NSK**

Please read this Operation Manual carefully before use, and file for future reference.

High Speed Air Turbine Handpiece

**S·Max pico**  
S-Max pico

**MADE IN JAPAN**

OM-T0497E

Powerful Partners®

## Intended to Use

This handpiece is designed only for dental treatment use.

### CAUTION

- When operating the Handpiece always consider the safety of the patient.
- Do not exceed maximum bur length.
- Check for vibration, noise and overheating outside the patient's oral cavity before use. If any abnormalities are found, stop using the Handpiece immediately and contact dealer.
- Should the handpiece function abnormally, cease operation immediately and return the handpiece to dealer for repair.
- Depressing the Push Button while the handpiece is rotating will result in OVERHEATING of the handpiece head. Special caution must be exercised during use to keep cheek tissue **AWAY** from the Push Button of the handpiece. Contact with cheek tissue may cause the Push Button to depress and burn injury to the patient may occur.
- Do not applying excessive power in the direction which bur removes. When the bur drops into patient's oral cavity, they may drink by mistake and get injured.
- Do not allow any impact on to the handpiece. Do not drop the handpiece.

- Do not attempt to disassemble the handpiece or tamper with the mechanism except NSK recommend on this Operation Manual.
- Do not connect or disconnect the handpiece or bur until the handpiece has completely stopped.
- NSK has not recommended the repair and disassembly of the cartridge without exception. Handpiece may perform abnormally (abnormal noise or abnormal vibration). Damage, failure or accidents are outside of our guarantee.
- Do not wipe with, or clean or immerse in, high acid water or sterilizing solutions.
- Use eye protectors, mask for safety when operating this handpiece.
- User is responsible for operation, maintenance and safety.
- If the handpiece has not been used for a long period, check for noise, vibration, and overheating before use.

# 1. Specification

Model	S-Max pico QD	S-Max pico M4	S-Max pico B2
Connection Type	QD Coupling	-	-
Hose Connection Type	-	ISO 9168 Type2 (Midwest 4 hole)	ISO 9168 Type1 (Borden 2 hole)
Rotation Speed	380,000 - 450,000 min <sup>-1</sup>		
Drive Air Pressure	0.20 MPa (2.0kgf/cm <sup>2</sup> ) - 0.25MPa (2.5kgf/cm <sup>2</sup> )	0.18 MPa (1.8kgf/cm <sup>2</sup> ) - 0.22MPa (2.2kgf/cm <sup>2</sup> )	
Spray Type	Single Spray		
Chuck Type	Push Button Chuck		
Bur	ISO1797-1 Type3 Ø1.59 -1.60mm Short Shank Bur / Super Short Shank Bur		
Max. Bur Length	18.5mm		
Bur Shank Length	8.0mm		
Max. Working Part Diameter	Ø 2.0mm		
Max.Air Consumption	40±5NL/min		
Head Size (external diameter x Hight)	Ø 8.6mm x 9.0mm		
Water Pressure	0.05MPa - 0.20MPa		
Chip Air Pressure	0.15MPa - 0.25MPa	-	
Use Environment	Temperature:10-40°C, Humidity: 30-75% RH, Atmospheric pressure: 700-1,060hPa		
Transportation and Storage Environment	Temperature:-10-50°C, Humidity: 10-85% RH, Atmospheric pressure: 500-1,060hPa		

## 2. Setting of Air & Water Supply Pressure

Measure the supply pressure at the handpiece / hose connection point and set the pressure to the value specified on the specification table. (Fig. 1)

\*Multi Gauge : Midwest 4 Hole Order code Z109400  
Borden2/3 Hole Order code Z109200


Fig. 1

### **WARNING**

Do not exceed the optimum pressure specified on the specification table.

### **CAUTION**

Do not use air contaminated by dust, moisture and oil.

### 3. Connecting/Disconnecting Handpiece

S-Max pico QD

#### (1) Connecting

- 1) Insert the Coupling into the Hose Connector.
- 2) Tighten the Hose Nut.
- 3) Pull the Connector Ring
- 4) Push into the handpiece into Coupling.

#### (2) Disconnect

Slide back the Connector Ring as shown in Fig. 2.


Fig. 2

#### CAUTION

- Ensure connect “Coupling and Hose” or “Handpiece and Coupling” firmly.
- Do not pull the Connector Ring while the handpiece is under air pressure. It may jump the handpiece out from the connector because of the air pressure.

## S-Max pico M4/S-Max pico B2

### (1) Connecting

- 1) Insert the handpiece correctly into the hose connector and tighten the Hose Nut. (Fig. 3)
- 2) Make sure the handpiece is connected firmly to the hose.

### (2) Disconnect

Loosen the Hose Nut and remove from the Hose.


Fig. 3

## 4. Insertion & Removal of the bur

### CAUTION

- Do not use any bur longer than recommended by NSK in Specification.
- Do not use bent, worn, damaged, or non-concentric burs. Such burs can cause damage to the handpiece.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Always keep the bur shank clean. Entry of hard debris into the chuck via the bur shank could cause rotation slip and also prevent the bur from being securely located in the chuck.

#### 4-1 To insert the Bur

- 1) Insert the bur until it stops. (①)
- 2) Depress the Push Button. (②)
- 3) Insert the bur into the chuck until it is stops. (③)
- 4) Release the button.

#### 4-2 To remove the Bur

- Depress the Push Button and remove the bur.  
(④)


Fig. 4

## CAUTION

- Fully depress the Push Button and insert the bur into the chuck until it is secure.
- Do not use carbide bur when the crown is removed. It may cause malfunction.
- Do not apply too much pressure to the bur. Bur may break or bent or difficult to remove. (Be sure to irrigate water. Cutting should be soft touch.)

## 5. Maintenances

### 5-1 Cleaning (Clean-Head System)

- 1) Half fill a cup with clean water.
- 2) Scrub dirt and debris from the handpiece, and wipe clean with alcohol-immersed cotton swab or cloth.  
(Fig. 5) Do not use a wire brush.
- 3) Rotate the handpiece and immerse half of the handpiece head. (Fig. 6)
- 4) Rotate and stop intermittently the handpiece about 2 to 3 seconds. (About 3 times)
- 5) By the function of the clean head system, the clean head hole can be washed.
- 6) Stop the handpiece and wipe it dry. If the dirt could not remove from the hole, remove the cartridge and clean it by brush.
- 7) Perform lubrication and sterilization. (5-4 Lubrication and 5-5 Sterilization)


Fig. 5


Fig. 6

### CAUTION

Use only clean water to maintain the Clean Head System.

### 5-2 Cleaning (Spray Port)

When Spray Port are clogged, clean the port as follows:

- 1) Remove the bur from the chuck.
- 2) Take out the cleaning wire from the back end of the brush holder. Insert the wire straight into the spray port with caution. Do not forcibly insert the wire into the port. Damaged or disoriented port could cause the spray diverge or directed away from the bur, and the cooling efficiency deteriorates.
- 3) After use, clean the cleaning wire. Push into the brush holder with the pointed end of the wire inward.


Fig. 7

### 5-3 Cleaning (Handpiece)

 This handpiece can be washed via Thermo Disinfector.

Scrub dirt and debris from the handpiece, and wipe clean with alcohol-immersed cotton swab or cloth. Do not use a wire brush.

\* Refer the manual of Thermo-Disinfector for performing it.

#### 5-4 Lubrication

##### ■ PANA SPRAY Plus /PANA SPRAY

Apply NSK PANA SPRAY Plus/PANA SPRAY/automatic handpiece cleaning and lubrication system every time after each use and/or before autoclaving.

##### S-Max pico QD

- 1) Remove the bur from the handpiece.
- 2) Insert the Spray Nozzle into the Spray Port nozzle on the can.
- 3) Insert the Spray Nozzle in rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds. (Fig. 8)

##### S-Max pico M4/S-Max pico B2

- 1) Remove the handpiece from the hose.
- 2) Mount the arrow-head spray nozzle tip into the spray can port. (Fig. 9)
- 3) Insert the Spray Nozzle into the Drive Air Port of the handpiece. (Fig. 10) Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds.


Fig. 8


Fig. 9

## CAUTION

- Be sure to hold the handpiece firmly to prevent the handpiece from slipping out of hand by the spray pressure when lubricating.
- Supply lubricant until it comes out of the handpiece head.(at least 2 seconds)
- Hold the spray can upright.

### ■ Lubrication Oil

S-Max pico M4/S-Max pico B2

- 1) Remove the handpiece from the hose.
- 2) Apply 1-2 drops of lubrication oil into the Drive Air Port shown in the Fig. 10.
- 3) Connect the handpiece to the hose. Operate the handpiece for about 10 seconds.


Fig. 10

### ■ Chuck cleaning

Clean Push Button chuck once a week.

- 1) Mount the arrow-head spray nozzle tip into the spray can port.
- 2) Lubricate the chuck directly through the bur insertion hole. (Fig. 11)
- 3) Lubricate the handpiece by using NSK PANA SPRAY Plus/PANA SPRAY (Fig. 8, 9) or NSK automatic handpiece cleaning and lubrication system.


Fig. 11

## CAUTION

If the chuck is not regularly cleaned the chuck grip may be weakened and the bur may be accidentally released while in operation.

### ■ NSK automatic handpiece cleaning and lubrication system

When using NSK automatic handpiece cleaning and lubrication system, refer to the system instructions.

## 5-5 Sterilization

Autoclave sterilization is recommended. Sterilization is required first time you use and after each patient as noted below.

### ■ Autoclave Procedure:

- 1) Scrub dirt and debris from the handpiece, and wipe clean with alcohol-immersed cotton swab or cloth. Do not use a wire rush.
- 2) Lubricate with PANA SPRAY Plus/PANA SPRAY/NSK automatic handpiece cleaning and lubrication system. Refer to 5-4 Lubrication.
- 3) Insert into an autoclave pouch. Seal the pouch.
- 4) Autoclavable up to max. 135°C.  
ex.) Autoclave for 20 min. at 121°C, or 15 min. at 132°C.
- 5) Keep the handpiece in the autoclave pouch to keep it clean until you use it.

\*Sterilization at 121°C for more than 15 minutes is recommended by EN13060 or EN ISO17665-1.

## CAUTION

- To prevent infectiveness, Autoclave the handpiece every patient.
- Do not autoclave the handpiece, even when it is bagged, together with other instruments, from which chemicals are not sufficiently cleaned. This is to avoid discoloration and/or damage to the plating of the handpiece.

- Do not wash, soak, or wipe off the handpiece with/in the oxidation potential solution (strong acid, super acid solution) sterilized solution, strong alkaline solution or functional water.
- The handpiece is kept in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. Moreover, air does not include a dust, salt and sulfur.
- Do not touch the handpiece just after the autoclave. It may be very hot.
- NSK do not recommend you to sterilize on plasma sterilization or EOG sterilization.

## 6. Replacing the Cartridge

- 1) Insert a dummy bur into the chuck.
- 2) Unscrew the Head Cap by turning counter-clockwise and remove.
- 3) Push up the dummy bur upwards to remove Cartridge.
- 4) Clean the interior of the Head casing.
- 5) Insert the new Cartridge.
- 6) Replace the Head Cap and firmly screw into place.

\*SX-PU03: Order Code P1140200


Fig. 12

## CAUTION

- Make sure that the O-Ring is removed, when removing the Cartridge. O-Ring might be remained in the Head, when removing the Cartridge. If there is a remained O-Ring, the new Cartridge can not be inserted properly.
- Use only NSK genuine cartridge. If you use other cartridge, specification may not be satisfied and we can not warranty of this situation.
- The head cap screw thread is very fine. To prevent damage to the thread always first tighten the handpiece head cap using fingers, then secure with the head cap wrench.
- Make sure that the O-ring (Fig. 12) is correctly located on the cartridge.

## 7. Replacing the O-rings (QD Coupling)

Worn or damaged O-ring may cause air or water to leak.

When water is found in the exhaust air or water leak from the coupled point, replace our O-rings as shown in Fig. 13.

\*Replacement O-ring Set (4pcs) : Order Code Y900592


Fig. 13 15

## 8. Non-Retraction Valve (QD Coupling)

To prevent oral fluids drawn into the water line, QD coupling has Non-Retraction valve. If water is beginning to leaking from handpiece, replace non-retraction valve.

- 1) Remove the back-end gasket.
- 2) Replace Non-Retraction Valve.


Fig. 13

## 9. Option Parts List

Model	PC1	PC2	PC3	PC4	PC5	PC6
Bur						
Max working length	16.5 mm	16.5 mm	17.5 mm	16.0 mm	14.0 mm	15.5 mm
Max working portion length	1.9 mm	0.4 mm	9.0 mm	7.0 mm	4.0 mm	6.0 mm
Working portion diameter	0.9 mm	0.6 mm	1.8 mm	1.3 mm	1.6 mm	1.1 mm
Grit	Standard	Fine	Coarse	Standard		
Order Code	Y1001817 (3pcs)	Y1001818 (3pcs)	Y1001819 (3pcs)	Y1001820 (3pcs)	Y1001821 (3pcs)	Y1001822 (3pcs)

Model	PC3ff	PC4ff	PC5ff	PC6ff
Bur				
Max working length	17.5 mm	16.0 mm	14.0 mm	15.5 mm
Max working portion length	9.0 mm	7.0 mm	4.0 mm	6.0 mm
Working portion diameter	1.7 mm	1.2 mm	1.6 mm	1.1 mm
Grit	Extra Fine			
Order Code	Y1002019 (3pcs)	Y1002020 (3pcs)	Y1002021 (3pcs)	Y1002022 (3pcs)

## 10. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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 The EU directive 93/42/EEC was applied in the design and production of this medical device.

---

**NAKANISHI INC.** 

[www.nsk-inc.com](http://www.nsk-inc.com)

700 Shimohinata Kanuma-shi  
Tochigi 322-8666,  
Japan

**NSK Europe GmbH** 

[www.nsk-europe.de](http://www.nsk-europe.de)

Elly-Beinhorn-Strasse 8  
65760 Eschborn,  
Germany

**NSK United Kingdom Ltd**

[www.nsk-uk.com](http://www.nsk-uk.com)

Office 5, Gateway1000,  
Arlington Business Park, Whittle Way,  
Stevenage, SG1 2FP, UK

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**NSK Oceania Pty Ltd**

[www.nskoceania.com.au](http://www.nskoceania.com.au)

Unit 22, 198-222 Young St.  
Waterloo, Sydney,  
NSW 2017, Australia

**NSK Middle East**

[www.nsk-inc.com](http://www.nsk-inc.com)

Room 6EA-701, 7th Floor, East Wing No.6  
Dubai Airport Free Zone,  
PO Box 54316 Dubai, UAE

**NSK Asia Pte Ltd**

[www.nsk-inc.com](http://www.nsk-inc.com)

1 Maritime Square,  
#09-33 HarbourFront Centre,  
Singapore 099253

Specifications are subject to change without notice.

2016.08.29 003<sup>®</sup>

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**NSK France SAS**