




NSK presents
**ultra mini head &
super slim body**

The NSK S-Max pico series has been development
for the operation of Minimal intervention (M.I).
NSK is ensure that wide visibilities and airy
operational feeling.



new

S-Max

High Speed Handpieces

pico



Powerful Partners®

S-Max pico


- Stainless Steel Body
- ISB Turbine
- Cellular Glass Optics
- Ceramic Bearings
- Clean Head System
- Single Spray

LED




NSK Couplings

Comparison of body sizes



Size specifications of the S-Max pico series



S-Max pico (Ultra Miniature Head)

Ti-Max X600L (Standard Head)



S-Max pico



Ti-Max X600L

Cartridge Replacement, Water Spray & Cellular Optics



Connect direct to major couplings

